



SP250

The ultimate way to convert your wheelchair to motorized wheelchair.

Allows you to keep the advantages of Manual-Foldable wheelchair as well as Powered wheelchair.

Use it for long trips, or remove the powered wheels and keep your good shape while using it as manual.

Fold it into a regular car trunk as powered or manual.

SPECIFICATION

| | |
|-------------------|--|
| Motor | 20/ 22/ 24-inch Brushless Gear motor + spoke |
| Motor weight | 7kg*2 |
| Brake | Electromagnetic brake |
| Speed | 1-5 level (Min. 3km/h, Max. 8km/h) |
| Rated power | 250W |
| Rated Voltage | 24V |
| No load current | < 0.6A |
| No load Speed | 100±10 rpm/min |
| Battery | 10AH |
| Operating current | 2-4A |
| Max. Torque | 70N.m |
| The Load Capacity | 120kg |
| Slope | ≤12° |

