

DRYFLOAT



DryFloAT Anti Decubitus non electric dynamic seating cushion.

DryFloAT offers an optimal anatomic immersion, low pressure, low humidity, low shear & friction as well as, high skin air ventilation environment, for individuals with skin/soft tissue integrity at all risk levels.

Pliable 5 to 10cm tall cells, allow optimal immersion into the cushion for maximum pressure redistribution, while decreasing friction and shear. The simple adjustment allows the cushions to fit the shape of the individual, minimizing tissue deformation and maximizing blood flow.



Whether you are looking for pain free seating, or had history of skin/soft tissue breakdown and/or any stage pressure injury (including deep tissue injury), you should try the DryFloAT Anti Decubitus non electric dynamic seating cushion

DryFloAT is adjustable, washable & repairable design and produced, based on 30 years of actual field experience as well as clinical studies, showing the medical advantages of the DryFloAT technology.

DryFloAT is available to fit individual measurements from 35cm to 55cm width & depth.

